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Referenzen:

1. Kornberg et al. (2016). Clinical Utility of Ultra-Widefield Imaging with the Optos Optomap Compared with Indirect Ophthalmoscopy in the Setting of Non-Traumatic Rhegmatogenous Retinal Detachment. *Seminars in ophthalmology*, 31(5), 505–512. <https://doi.org/10.3109/08820538.2014.981551>
2. Abadia et al. (2020). Non-Mydriatic Ultra-Wide Field Imaging Versus Dilated Fundus Exam and Intraoperative Findings for Assessment of Rhegmatogenous Retinal Detachment. *Brain sciences*, 10(8), 521. <https://doi.org/10.3390/brainsci10080521>
3. Bacherini et al. (2025). Characterization of Peripheral Retinal Degenerations and Rhegmatogenous Lesions Using Ultra-Widefield Swept Source OCT Integrated with a Novel Scanning Laser Ophthalmoscope. *Diagnostics* (Basel, Switzerland), 15(22), 2930. <https://doi.org/10.3390/diagnostics15222930>
4. Lin et al. (2021). Ultra-widefield retinal imaging for adjunctive resident training in retinal break detection. *PLoS one*, 16(6), e0253227. <https://doi.org/10.1371/journal.pone.0253227>
5. McGrath et al. (2025). Sensitivity and utility of ultra-wide field imaging for the detection of peripheral retinal breaks in a large Irish tertiary referral centre. *Eye* (London, England), 39(16), 2980–2985. <https://doi.org/10.1038/s41433-025-03998-6>
6. Francone et al. (2020). DETECTION OF NEUROSENSORY RETINAL DETACHMENT COMPLICATING DEGENERATIVE RETINOSCHISIS BY ULTRA-WIDEFIELD FUNDUS AUTOFLUORESCENCE IMAGING. *Retina* (Philadelphia, Pa.), 40(5), 819–824. <https://doi.org/10.1097/IAE.0000000000002488>
7. Zhang et al. (2025). Ultra-Widefield Fundus Autofluorescence Findings as an Indicator of Duration in Rhegmatogenous Retinal Detachment. *Ophthalmology science*, 6(2), 101016. <https://doi.org/10.1016/j.xops.2025.101016>
8. Poulsen et al. (2019). Fundus autofluorescence and spectral domain optical coherence tomography as predictors for long-term functional outcome in rhegmatogenous retinal detachment. *Graefes archive for clinical and experimental ophthalmology = Albrecht von Graefes Archiv für klinische und experimentelle Ophthalmologie*, 257(4), 715–723. <https://doi.org/10.1007/s00417-018-04222-w>
9. Ishikawa et al. (2019). Preoperative estimation of distance between retinal break and limbus with wide-field fundus imaging: Potential clinical utility for conventional scleral buckling. *PLoS one*, 14(2), e0212284. <https://doi.org/10.1371/journal.pone.0212284>
10. Witmer et al. (2012). Ultra-wide-field autofluorescence imaging in non-traumatic rhegmatogenous retinal detachment. *Eye* (London, England), 26(9), 1209–1216. <https://doi.org/10.1038/eye.2012.122>
11. Lasave et al. (2019). When Ultra-Widefield Retinal Imaging Is Most Useful. *Retina Today*, Issue April 2019 Pages 36-38, 40, 12. Li et al. (2023). Preoperative estimation of retinal hole location using ultra-wide-field imaging. *Annals of medicine*, 55(2), 2258790. <https://doi.org/10.1080/07853890.2023.2258790>
13. Mahmoudzadeh et al. (2024). Retinal Displacement Following Rhegmatogenous Retinal Detachment Repair. *Ophthalmic surgery, lasers & imaging retina*, 55(10), 560–566. <https://doi.org/10.3928/23258160-20240528-01>
14. Valentin-Bravo et al. (2024). Silicone oil emulsification: A literature review and role of widefield imaging and ultra-widefield imaging with navigated central and peripheral optical coherence tomography technology. *Saudi journal of ophthalmology : official journal of the Saudi Ophthalmological Society*, 38(2), 112–122. https://doi.org/10.4103/sjopt.sjopt_193_23

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